

Princeton Chromatography Inc.

Princeton Chromatography

PrincetonSFC 2CN-DIOL



PrincetonSFC-2CN/DIOL 60Å..., 100Å... pore size and 5 µm particles Bonded Phase HPLC columns are prepared with a high quality spherical silica.

The low surface areas of these packing materials are beneficial for separations requiring low percentages of organic modifiers.

The carbon load for this material is 16%(60Å...) and 9%(100Å...)



075046-03086-S
\$270.00

[\[Product Details...\]](#)



075046-03586-S
\$345.00

[\[Product Details...\]](#)



075070-01586-S
\$375.00

[\[Product Details...\]](#)



075070-03586-S
\$375.00

[\[Product Details...\]](#)



100030-01586-S
\$370.00

[\[Product Details...\]](#)



100030-03586-S
\$370.00

[\[Product Details...\]](#)



100040-01086-S
\$310.00

[\[Product Details...\]](#)



100040-01586-S
\$370.00

[\[Product Details...\]](#)



100040-03086-S
\$280.00

[\[Product Details...\]](#)



100040-03586-S
\$370.00

[\[Product Details...\]](#)



100046-01086-S
\$310.00

[\[Product Details...\]](#)



100046-01586-S
\$370.00

[\[Product Details...\]](#)



100046-03086-S
\$280.00

[\[Product Details...\]](#)



100046-03586-S
\$370.00

[\[Product Details...\]](#)



100070-01586-S
\$400.00

[\[Product Details...\]](#)



100070-03586-S
\$400.00

[\[Product Details...\]](#)



100078-01086-S
\$625.00

[\[Product Details...\]](#)



100078-01586-S
\$735.00

[\[Product Details...\]](#)



100078-03086-S
\$625.00

[\[Product Details...\]](#)



100078-03586-S
\$735.00

[\[Product Details...\]](#)

-
-
-
-
-
-
-
-
-
-
-

[« « Start](#)

[« Prev](#)

[1](#)

[2](#)

[3](#)

[4](#)

[5](#)

[6](#)

[7](#)

[Next »](#)

[End »»](#)

Results 41 - 60 of 140

Recently Viewed Products

- [050078-01079-S](#) (Category: [PrincetonSemi-Prep SFC](#))
- [050078-01586-S](#) (Category: [PrincetonSFC 2CN-DIOL](#))
- [050078-03086-S](#) (Category: [PrincetonSFC 2CN-DIOL](#))
- [050030-01551-H](#) (Category: [HPLC Columns](#))
- [050030-03309-S](#) (Category: [PrincetonSFC Analytical](#))