



**Column:** PrincetonSPHER 2-Ethylpyridine  
**Catalog No:** 150046-01577  
**Dimensions:** 150 x 4.6 mm  
**Eluent:** 80/20 A/B

A: 20mM K<sub>2</sub>HPO<sub>4</sub> - pH 6.6  
B: Acetonitrile

**Flow Rate:** 1.0 mL/min  
**Temperature:** Ambient  
**Detection:** UV at 273 nm  
**Injection:** 10uL (1 mg/mL sample)

1. Metoprolol
2. Propanolol
3. Pindolol