



Column: PrincetonSPHER C18 100A 5u
Catalog No: 150046-03501
Dimensions: 150 x 4.6 mm
Eluent: 60/40 A/B

A: 20mM NH₂H₂PO₄ - pH 3.2
B: Acetonitrile

Flow Rate: 1.5 mL/min
Temperature: 40 degrees C
Detection: UV at 280 nm
Injection: 10 uL (1 mg/mL sample)

1. Hydrochlorothiazide
2. Chlorthalidone
3. Diazoxide
4. Furosemide
5. Bendroflumethazide