



**Column:** PrincetonSPHER OmniSEP  
**Catalog No:** 250046-01590  
**Dimensions:** 250 x 4.6 mm  
**Eluent:** 0.1% H<sub>3</sub>PO<sub>4</sub>

**Flow Rate:** 1.0 mL/min  
**Temperature:** 30 degrees C  
**Detection:** UV at 210 nm  
**Injection:** 20 uL (1 mg/mL sample)

1. Tartaric Acid
2. Formic Acid
3. Malic Acid
4. Malonic Acid
5. Lactic Acid
6. Maleic Acid