



Column: PrincetonSPHER Ultima-C18-WP
Catalog No: 150046-01524
Dimensions: 150 x 4.6 mm
Eluent: 40/60 A/B

A: 25mM KH₂PO₄ - pH 6.6
B: Acetonitrile

Flow Rate: 2.0 mL/min
Temperature: 40 degrees C
Detection: UV at 230 nm
Injection: 10uL (1 mg/mL sample)

1. Protriptyline
2. Desipramine
3. Nortriptyline
4. Doxepin
5. Imipramine
6. Amitriptyline
7. Trimipramine