



**Column:** PrincetonSPHER C18 100A 5u  
**Catalog No:** 100046-03501  
**Dimensions:** 100 x 4.6 mm  
**Eluent:** 50/50 A/B

A: 0.1M H<sub>3</sub>PO<sub>4</sub> - pH 2.  
B: Methanol

**Flow Rate:** 1.5 mL/min  
**Temperature:** Ambient  
**Detection:** UV at 240 nm  
**Injection:** 5 uL (1 mg/mL sample)

1. p-Aminobenzoic Acid
2. p-Hydroxybenzoic Acid
3. o-Aminobenzoic Acid
4. Benzyl Alcohol
5. Salicylic Acid