



**Column:** PharmaBOND C18 125A 10u  
**Catalog No:** 300039-05001  
**Dimensions:** 300 x 3.9 mm  
**Eluent:** 60/36/4; Methanol, water, Tetrahydrofuran w/ 0.1% H3PO4 and 4.4mg/mL Doc

**Flow Rate:** 1.5 mL/min  
**Temperature:** Ambient  
**Detection:** UV at 280 nm  
**Injection:** 15uL (1 mg/mL sample)  
**Sample:** 0.85 g Acetaminophen/ Codeine tablet (#3) in 250 mL mobile phase

1. Acetaminophen
2. Codeine