



Column: PrincetonSPHER FO* 60A 5u
Catalog No: 150046-01542
Dimensions: 150 x 4.6 mm
Eluent: 30/70 A/B

A: Acetonitrile
B: 0.1% H₃PO₄

Flow Rate: 1.0 mL/min
Temperature: 20degrees C
Detection: UV at 254 nm
Injection: 5 uL (1 mg/mL sample)

1. Vanillic Acid
2. Vanillin
3. Ethyl Vanillin

*FO = Fluorooctyl