



Column: PrincetonSPHER C18 60A 5u
Catalog No: 250046-01501
Dimensions: 250 x 4.6 mm
Eluent: 70/30 A/B

A: Acetonitrile
B: Water

Flow Rate: 1.5 mL/min
Temperature: 30 degrees C
Detection: UV at 254 nm
Injection: 10 uL (1 mg/mL sample)

1. 2,4-DNT
2. 4,4'-DPA
3. 4N-DPA
4. NNO-DPA
5. 2,4-DPA
6. 2-2'-DPA
7. 2,4-DPA
8. DPA
9. 2N-DPA
10. Dibutyl Phthalate